



Philips
USB Flash Drive

8GB
blue edition

FM08FD02B



Share and transport data

Swap files in a flash!

The Philips USB flash drive blue edition is an easy way to share and transport your data. A useful 8GB gives you the capacity to swap and share large files like documents, pictures and music via the USB port of your PC or laptop.

Everything you need

- Specially designed and easy to use software pre-installed
- Activity indicator lights up when copying files

Performance and capacity

- 8GB storage capacity for large data files
- Automatic compression of all files saves free space
- Fast data transfer with high-speed USB 2.0
- File synchronisation between complete PC and UFD

Safe and secure

- Lost & Found details so the drive can be returned if lost
- Protect your valuable, personal data with a password

PHILIPS
sense and simplicity

USB Flash Drive
8GB blue edition

Specifications

Storage Media

- **Built-in memory type:** MLC NAND Flash
- **Built-in Memory Capacity:** 8 GB
- **Transfer rate:** reads max 15MB/s, write max 3MB/s

Connectivity

- **USB:** High-speed USB 2.0

Software

- **Carry it Easy:** Password protection, File synchronisation, Portable Outlook Express-trial, Favourites Sync-trial, Trace free browsing-trial

Accessories

- **Neck strap:** No
- **Cables:** None
- **Quick start guide:** English, French, German, Spanish, Italian, Dutch, Portuguese, Russian
- **User Manual:** English, French, German, Spanish, Italian, Dutch, Portuguese, Russian

System Requirements

- **PC OS:** Windows® Vista, XP, 2000, ME and 98SE; Mac OS 9.0 and higher; Linux 2.4.0 and higher
- **USB:** Free USB port

Miscellaneous

- **Warranty:** 2 year

Packaging Data

- **I2NC:** 8670 000 32851
- **EAN/UPC/GTIN:** 87 12581 35572 2
- **Height:** 180 mm
- **Length:** 130 mm
- **Width:** 15 mm
- **Gross weight:** 0.048 kg
- **Tare weight:** 0.039 kg

Inner Carton

- **Quantity:** 15
- **EAN/UPC/GTIN:** 87 12581 35573 9
- **Length:** 200 mm
- **Height:** 140 mm
- **Gross weight:** 0.86 kg
- **Tare weight:** 0.725 kg
- **Width:** 180 mm

Outer Carton

- **Quantity:** 60
- **EAN/UPC/GTIN:** 87 12581 35574 6
- **Length:** 315 mm
- **Width:** 230 mm
- **Height:** 395 mm
- **Gross weight:** 4.052 kg
- **Tare weight:** 3.512 kg

FM08FD02B/00

Product highlights

8GB storage capacity

A useful 8GB gives you the capacity to swap and share much larger or more files via the USB port of your PC or laptop.

Automatic file compression

Save free space by automatic compression of all the files you put on the USB flash drive

Fast data transfer

High-speed data transfer rates significantly cut annoying waiting time when you're copying large multimedia files to or from your computer's hard disk.

File synchronisation

With file synchronisation of all the files on your PC, it is possible to access and work with your documents anywhere with just a few mouse clicks.

Lost & Found information

This unique feature enables you to enter your contact details like address, telephone number and e-mail address on the USB flash drive so that if you lose it, whoever finds it can send it back to you.

Password protection

Due to password protection including 256-bit AES encryption, your valuable personal data can be protected from unauthorised access by simply using a password.

Activity indicator

An appealing glow shows you the drive is properly connected and working and pulsates faster when copying files to or from its memory.



Issue date 2007-12-14

Version: 2.0.4

I2 NC: 8670 000 32851
EAN: 87 12581 35572 2

© 2007 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com