

Philips
USB Flash Drive

4GB

vintage edition

FM04FD20B



Designed to share, store and protect

your valuable data

Good looks, high performance and thoughtful design in the palm of your hands. The Philips USB flash drive vintage edition features a sleek design and useful touches like a captive cap adding extra appeal.

Everything you need

- Specially designed and easy to use software pre-installed

Performance and capacity

- Fast data transfer with high-speed USB 2.0
- Protect your valuable, personal data with a password
- 4GB storage capacity for large data files

Designed for you

- Activity indicator lights up when copying files
- Trendy, exclusive design with nice shape and smooth curves
- Integrated cap stays attached to the drive, never lose it!

PHILIPS
sense and simplicity

USB Flash Drive
4GB vintage edition

Specifications

Storage Media

- **Built-in Memory Capacity:** 4 GB
- **Transfer rate:** reads max 15MB/s, writes max 3MB/s
- **Built-in memory type:** MLC NAND Flash

Connectivity

- **USB:** High-speed USB 2.0

Software

- **Carry it Easy:** Password protection, File synchronisation, Portable Outlook Express-trial, Trace free browsing-trial, Favourites Sync-trial

Accessories

- **Cables:** No
- **Neck strap:** No
- **Quick start guide:** English, French, German, Spanish, Italian, Dutch, Portuguese, Russian
- **User Manual:** English, French, German, Spanish, Italian, Dutch, Portuguese, Russian

System Requirements

- **PC OS:** Windows® Vista, XP, 2000, ME, 98 SE
Mac OS 9.0 and higher; Linux 2.4.0 and higher
- **USB:** Free USB port

Miscellaneous

- **Warranty:** 2 year

Packaging Data

- **I2NC:** 8670 000 31612
- **EAN/UPC/GTIN:** 87 12581 34677 5
- **Height:** 180 mm
- **Length:** 130 mm
- **Width:** 15 mm
- **Gross weight:** 0.053 kg
- **Tare weight:** 0.039 kg

Inner Carton

- **Quantity:** 15
- **EAN/UPC/GTIN:** 87 12581 35580 7
- **Length:** 200 mm
- **Width:** 180 mm
- **Height:** 140 mm
- **Gross weight:** 0.935 kg
- **Tare weight:** 0.725 kg

Outer Carton

- **EAN/UPC/GTIN:** 87 12581 35581 4
- **Length:** 315 mm
- **Width:** 230 mm
- **Height:** 395 mm
- **Gross weight:** 4.352 kg
- **Tare weight:** 3.512 kg

FM04FD20B/00

Product highlights

Activity indicator

An appealing glow shows you the drive is properly connected and working and pulsates faster when copying files to or from its memory.

Exclusive design

With its smooth curves and a soft-touch rubber finish, the Philips USB flash drive is an exclusive device with a stunning design.

Fast data transfer

High-speed data transfer rates significantly cut annoying waiting time when you're copying large multimedia files to or from your computer's hard disk.

Password protection

Due to password protection including 256-bit AES encryption, your valuable personal data can be protected from unauthorised access by simply using a password.

Easy to use software included

A complete software suite is included, which will run automatically from the USB flash drive and gives you access to useful features like, Password protection, File synchronisation, Data compression, USB lost & found and USB information. Additionally, it comes with a 30-day trial for Portable Microsoft Outlook Express, Favourites synchronisation and No-trace internet browsing.

Integrated cap

The cap to protect the USB connector is integrated in the product design so that it stays attached to the drive when not in use and is impossible to lose.

4GB storage capacity

A useful 4GB gives you the capacity to swap and share much larger or more files via the USB port of your PC or laptop.



Issue date 2007-12-14

Version: 2.0.4

I2 NC: 8670 000 31612
EAN: 87 12581 34677 5

© 2007 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com